SURFACES & EQUIPMENT ANALYSIS SWABS TAMPONI PER ANALISI SUPERFICI

Swabs with liquid neutralizing solution, designed for microbiological surface analysis in food and pharmaceutical industries. Swabs with a plastic shaft with break point and a tip in polyurethane foam or viscose, for easier surface sampling. Tube Ø16x100 mm, with support base and screw cap "capture" type, containing a well-defined amount of specific solution to neutralize the possible presence of disinfectants or sanitizers on surfaces. Devices in single peel-pack, sterile (SAL 10⁻⁶), suitable for sampling in compliance with international standards ISO 18593, ISO 17604, ISO 6579, ISO 11290-1, ISO 11290-2, ISO 21528-1 and ISO 6887.

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SWAB + TEST TUBE

TAMPONE + PROVETTA


Swab with plastic shaft with break-point and Ø 16 x 100 mm conical tube, equipped with a self-standing base and yellow screw cap, containing 10 ml or 5 ml of neutralizing saline inhibitory isotonic solution. The tube is pre-labeled to ensure proper identification of the sample and collection site.

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SWAB PRE-INSERTED IN TEST TUBE

TAMPONE PRE-INSERITO IN PROVETTA


Swab with plastic shaft pre-inserted in the tube, already attached to the yellow screw cap. Ø 16 x 100 mm conical tube with self-standing base, containing 10 ml or 5 ml of neutralizing saline inhibitory isotonic solution. The tube is pre-labeled to ensure proper identification of the sample and collection site.

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Sampling surface area template made of honeycomb polypropylene, thick 2.5 mm. Available in versions with a sampling area of 10 x 10 cm (100 cm²) or 5 x 4 cm (20 cm²). Sterilized by irradiation.

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SURFACES AND EQUIPMENT ANALYSIS SPONGE

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NEW
 product

Sterile sponges for microbiological sampling according to European regulations and international standards ISO 18593, ISO 17604, ISO 6579, ISO 11290, ISO 21528-1 and ISO 6887, specifically intended for the analysis of surfaces, equipment, and animal carcasses to detect the presence of *Listeria* spp., *Staphylococci*, *Enterobacteriaceae*, *Salmonella*, *E. coli*, *Enterobacter*, *Klebsiella*, and other pathogens. The medical-grade polyurethane foam sponge, measuring 4 x 7 cm (1.5 cm thick), does not leave residues, does not break, and does not drip. The sponge is contained in a 16 x 22 cm nylon (PA/PE) bag with a strong thickness (70 microns), equipped with a tear seal (to open it, it is necessary to tear along the dotted line at the top to ensure sterility) and a zip-lock closure system that is liquid-tight. Each bag shows proper instructions use and a writing area for suitable sample identification. The bag can be used for pre-enrichment in the search of pathogenic bacteria by direct adding of appropriate selective culture broth. Additionally, it can be turned inside out and used as a glove during sampling and can be used in a homogenizer.

Sterile sponges can be supplied in kits in combination with single-use sterile sampling surface template measuring 10 x 10 cm, sterile gloves, and different solutions according to the sampling to be performed:

- **NEUTRALIZING DILUENT** (ISO 18593), the diluent solution with neutralizers inhibits the bacteriostatic properties of samples to undergo microbiological analysis.
- **MAXIMUM RECOVERY DILUENT** (Peptone Saline Water) is a protective isotonic diluent for optimal recovery of microorganisms from diversified sources (suitable for surfaces, foods, and carcasses) for qualitative and quantitative research (ISO 18593, ISO 6887). No neutralizing agents are included. For areas where disinfectant residues are suspected, the use of an alternative medium is recommended.
- **BUFFERED PEPTONE WATER**, a specific solution for animal carcasses (ISO 17604, ISO 6887, ISO 6579, ISO 11290, ISO 21528-1). It is recommended for the detection and count of *Salmonella* spp., *Enterobacteriaceae*, and *Listeria* spp. in foods. No neutralizing agents are included. For areas where disinfectant residues are suspected, the use of an alternative medium is recommended.

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
2472/SG/TL/G

2473/SG/TL/G

Swabs with selective culture medium for the identification of bacterial contamination on surfaces in food and pharmaceutical industries. Swabs with a plastic shaft and rayon tip, equipped with a cylindrical tube Ø 12 x 150 mm containing a chromogenic selective medium. After sampling, the swab is immersed in the tube's medium, and after an incubation period of 18/24 hours, the color change indicates the presence/absence of the contamination being tested. The label on the tube has a color legend for correct result interpretation. Devices in single peel-pack, sterile (SAL 10⁻⁶).

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XLD AGAR

Swabs with test tube containing XLD Agar (Xylose lysine deoxycholate agar), selective growth medium used for the isolation of Salmonella, Shigella and Coliform bacteria identification.

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2310/SG

OXFORD AGAR

Swabs with test tube containing OXFORD Agar, selective and diagnostic medium for the detection of Listeria monocytogenes.

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2320/SG

PLATE COUNT AGAR

Swabs with test tube containing PLATE COUNT Agar (Tryptic glucose yeast agar), suitable for the detection of total microbial contamination (aerobic and facultative anaerobes bacteria, heterophiles), for the bacteriological test on food, water, milk and other dairy products.

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2330/SG

COLIFORM AGAR

Swabs with test tube containing COLIFORM Agar, selective chromogenic medium for the detection of Escherichia coli and coliforms in food and water.

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