



CARRELLO INOX
Carrelli in acciaio inox AISI 304 a piani stampati. Struttura in tubo tondo di acciaio inox AISI 304 Ø cm 2,5, finitura spazzolata. Piani stampati a forma di vassoio con superfici arrotondate, insonorizzati, imbullonati al telaio. Paracolpi in gomma grigia. Portata kg 120. Ruote piroettanti Ø cm 12,5 di cui 2 dotate di freno. Viene fornito smontato.

ST. STEEL CART
Entirely made of AISI 304 stainless steel, brushed finish. Round tube frame, Ø 2.5cm. Tray shaped stainless steel pressed shelves bolted to frame, sound-proof and with no sharp edges. Grey rubber ring buffers. Capacity: 120 Kg. Swivel castors, Ø 12.5cm, 2 with brake. It is supplied disassembled.

Cod.	Lungh.-largh. Lenght-Width	Lungh.-largh.-alt. Leng.-Wid.-Heig.	Modello Model
05.1910.00	Vassoio\Trays 800x500	890x590x950	2 Vassoio\Trays
05.1911.00	Vassoio\Trays 800x500	890x590x950	3 Vassoio\Trays
05.1912.00	Vassoio\Trays 800x500	890x590x1260	4 Vassoio\Trays
05.1913.00	Vassoio\Trays 800x500	890x590x1600	5 Vassoio\Trays
05.1914.00	Vassoio\Trays 1000x500	1090x590x950	2 Vassoio\Trays
05.1915.00	Vassoio\Trays 1000x500	1090x590x950	3 Vassoio\Trays
05.1916.00	Vassoio\Trays 1000x500	1090x590x950	4 Vassoio\Trays
05.1917.00	Vassoio\Trays 1000x500	1090x590x950	5 Vassoio\Trays
05.1918.00	Vassoio\Trays 1000x600	1090x690x950	2 Vassoio\Trays
05.1919.00	Vassoio\Trays 1000x600	1090x690x950	3 Vassoio\Trays
05.1920.00	Vassoio\Trays 1000x600	1090x690x950	4 Vassoio\Trays
05.1921.00	Vassoio\Trays 1000x600	1090x690x950	5 Vassoio\Trays
05.1922.00	Vassoio\Trays 1200x600	1290x690x950	2 Vassoio\Trays
05.1923.00	Vassoio\Trays 1200x600	1290x690x950	3 Vassoio\Trays
05.1924.00	Vassoio\Trays 1200x600	1290x690x950	4 Vassoio\Trays



CARRELLO INOX CON PIANI ASPORTABILI
Carrelli in acciaio inox con piani asportabili. Struttura in tubo tondo di acciaio inox AISI 304 Ø cm 2,5, finitura spazzolata. Piani in lamiera di acciaio inox AISI 304, a forma di vassoio, con superfici arrotondate, profondi cm 3, estraibili. Paracolpi in gomma grigia. Portata kg 120. Ruote piroettanti Ø cm 12,5, di cui 2 dotate di freno. Viene fornito smontato.

ST. STEEL CART WITH REMOVABLE TRAY
AISI 304 stainless steel trolleys, brushed finish. Round tube frame Ø 2.5 cm. Removable soundproof tray-shaped shelves, with no sharp edges. Grey rubber ring buffers. Capacity: 120 Kg. Swivel castors, Ø 12.5 cm, 2 with brake. It is supplied disassembled.

Cod.	Lungh.-largh. Lenght-Width	Lungh.-largh.-alt. Leng.-Wid.-Heig.	Modello Model
05.1925.00	Vassoio\Trays 700x500	980x620x940	2 Vassoio\Trays
05.1926.00	Vassoio\Trays 700x500	980x620x940	3 Vassoio\Trays
05.1927.00	Vassoio\Trays 900x500	1180x620x940	2 Vassoio\Trays
05.1928.00	Vassoio\Trays 900x500	1180x620x940	3 Vassoio\Trays