


05.1170.00 - 05.1189.00



PORTAPROVETTE IN POLIPROPILENE AUTOCLAVABILI

Possono essere montati e smontati facilmente.

POLYPROPYLENE TEST TUBE RACK - AUTOCLAVABLE

This stand can be easily disassembled and reassembled.

Cod.	Diametro Diameter	N. buchi N. of holes	N. file Line n.	Colore Color
05.1170.00	13 mm	90	6x15	Bianco\White
05.1171.00	16 mm	60	5x12	
05.1172.00	20 mm	40	4x10	
05.1173.00	25 mm	24	3x8	
05.1174.00	30 mm	21	3x7	
05.1175.00	13 mm	90	6x15	Rosso\Red
05.1176.00	16 mm	60	5x12	
05.1177.00	20 mm	40	4x10	
05.1178.00	25 mm	24	3x8	
05.1179.00	30 mm	21	3x7	
05.1180.00	13 mm	90	6x15	Blu\Blue
05.1181.00	16 mm	60	5x12	
05.1182.00	20 mm	40	4x10	
05.1183.00	25 mm	24	3x8	
05.1184.00	30 mm	21	3x7	
05.1185.00	13 mm	90	6x15	Giallo\Yellow
05.1186.00	16 mm	60	5x12	
05.1187.00	20 mm	40	4x10	
05.1188.00	25 mm	24	3x8	
05.1189.00	30 mm	21	3x7	

05.1200.00 - 05.1203.00



PORTAPROVETTE IN POLIPROPILENE A TRE PIANI - AUTOCLAVABILI


Portaprovette in polipropilene a tre piani - autoclavabili.

THREE-LEVEL POLYPROPYLENE TEST TUBE RACK WITH SIDE HANDLES

Three-level polypropylene test tube rack with side handles.

Cod.	Diametro Diameter	Lungh.-largh.-alt. Leng.-Wid.-Heig.	N. buchi N. of holes	N. file Line n.
05.1200.00	13 mm	65x260x105 H. mm	31	3x10
05.1201.00	16 mm			
05.1202.00	20 mm	70x300x125 H. mm	20	2x10
05.1203.00	25 mm			

05.1250.00 - 05.1253.00



PORTAPROVETTE A 3 PIANI CON MANIGLIE LATERALI

In polipropilene - Autoclavabili. Resistenti alle basse temperature -80 °C / +125 °C.

3-LEVEL TEST TUBE RACKS WITH SIDE HANDLES

Made of polypropylene - Autoclavable. Resistant to low temperatures -80 °C / +125 °C.

Cod.	Diam. foro Hole diam.	Lungh.-largh.-alt. Leng.-Wid.-Heig.	N. buchi N. of holes	N. file Line n.
05.1250.00	13 mm	65x260x105 mm	31	3x10
05.1251.00	16 mm			
05.1252.00	20 mm	70x300x125 mm	20	2x10
05.1253.00	25 mm			

05.1350.00 - 05.1361.00



PORTAPROVETTE IN ACCIAIO INOX A 3 PIANI


Portaprovette in acciaio inox a 3 piani.

3-LEVEL STAINLESS STEEL TEST TUBE RACKS

3-Level stainless steel test tube racks.

Cod.	Diam. foro Hole diam.	Lungh.-largh.-alt. Leng.-Wid.-Heig.	N. buchi N. of holes	N. file Line n.
05.1350.00	13 mm	50x143x105 mm	12	2x6
05.1351.00	20 mm	65x175x105 mm		
05.1352.00	13 mm	74x143x105 mm	18	3x6
05.1353.00	20 mm	90x175x105 mm		
05.1354.00	13 mm	74x190x105 mm	24	3x8
05.1355.00	20 mm	90x230x105 mm		
05.1356.00	13 mm	120x143x105 mm	30	5x6
05.1357.00	20 mm	147x175x105 mm		
05.1358.00	13 mm	97x210x105 mm	36	4x9
05.1359.00	20 mm	120x260x105 mm		
05.1360.00	13 mm	97x280x105 mm	48	4x12
05.1361.00	20 mm	129x343x205 mm		

05.1800.00 - 05.1809.00



PORTAPROVETTE IN FILO DI ACCIAIO INOX


Portaprovette a 3 piani in filo inox

STAINLESS STEEL TEST TUBE RACK

Test tube rack at 3 shelves - PVC covered

Cod.	Diam. foro Hole diam.	Lungh.-largh.-alt. Leng.-Wid.-Heig.	N. buchi N. of holes	N. file Line n.
05.1800.00	13 mm	65x50x50 mm	12	4x3
05.1801.00	16 mm	120x45x65 mm		6x2
05.1802.00	25 mm	180x65x95 mm		
05.1803.00	16 mm	160x65x65 mm	24	8x3
05.1804.00	25 mm	235x95x95 mm		
05.1805.00	16 mm	160x125x65 mm	48	8x6
05.1807.00	25 mm	235x180x95 mm		
05.1808.00	13 mm	162X84X65 mm	50	10x5
05.1809.00	18 mm	210x110x65 mm		

05.1800.01 - 05.1807.01



PORTAPROVETTE RIVESTITO IN PVC

Portaprovette a 3 piani in filo inox rivestito in PVC

TEST TUBE RACK PVC COVERED

Test tube rack at 3 shelves - made of stainless steel - PVC covered

Cod.	Diam. foro Hole diam.	Lungh.-largh.-alt. Leng.-Wid.-Heig.	N. buchi N. of holes	N. file Line n.
05.1800.01	13 mm	65x50x50 mm	12	4x3
05.1801.01	16 mm	120x45x65 mm		6x2
05.1802.01	25 mm	180x65x95 mm		
05.1803.01	16 mm	160x65x65 mm	24	8x3
05.1804.01	25 mm	235x95x95 mm		
05.1805.01	16 mm	160x125x65 mm	48	8x6
05.1807.01	25 mm	235x180x95 mm		

05.1363.00 - 05.1366.00



PORTAPROVETTE CIRCOLARE PER BICCHIERI

Disponibile sia in acciaio inox che in alluminio.

ROUND TEST TUBE RACKS FOR BEAKERS

Available in both stainless steel and aluminum.

Cod.	Materiale Material	Diam. foro Hole diam.	Diametro 1 Diameter 1	Diametro 2 Diameter 2	N. buchi N. of holes
05.1363.00	Alluminio\ Aluminum	13-18 mm	Ø sup. 98,5 mm	Ø inf. 73,5 mm	4+4
05.1364.00	Acc. inox\ St. steel				
05.1365.00	Alluminio\ Aluminum	18 mm	Ø sup. 73,5 mm	Ø inf. 73,5 mm	8
05.1366.00	Acc. inox\ St. steel				

05.1368.00




PORTAPROVETTE CIRCOLARE IN POLIPROPILENE

In polipropilene - Autoclavabile. La base ha pozzetti semisferici coordinati e 16 perni verticali per asciugare i tubi circondati da un trogolo per raccogliere i gocciolamenti. Questo supporto può essere facilmente smontato e riasssemblato. Dispone di 4 posti con Ø 25 mm, 8 posti con Ø 19 mm e 16 supporti a pioli per scolo delle provette. Base stabile Ø 190 mm, con raccolta liquidi. Altezza totale 90 mm.

POLYPROPYLENE ROUND TEST TUBE RACK

Made of polypropylene - Autoclavable. The base has matching hemispherical wells and 16 verticals pins for drying tubes surrounded by a trough to catch drippings. This stand can be easily disassembled and reassembled. It has 4 places with Ø 25 mm, 8 places with Ø 19 mm and 16 upside down peg tube supports. Stable base Ø 190 mm, with liquid collection. Total height 90 mm.

05.1190.00 - 05.1191.00



PORTAPROVETTE IN POLIPROPILENE

In polipropilene - Autoclavabile, dotata di pioli per scolo delle provette. Questo supporto può essere facilmente smontato e riasssemblato.

POLYPROPYLENE TEST TUBE RACK

Made of polypropylene - Autoclavable, equipped with upside down peg tube supports. This stand can be easily disassembled & reassembled.

Cod.	Diametro 1 Diameter 1	Diametro 2 Diameter 2	Lungh.-largh.-alt. Leng.-Wid.-Heig.
05.1190.00	2 fori\2 holes Ø 25 mm	4 fori\4 holes Ø 16 mm	70x200x80 mm
05.1191.00	6 fori\6 holes Ø 25 mm	-	

05.1195.00



PORTAPIPETTE E TERMOMETRI ORIZZONTALE

In polipropilene - Autoclavabile. 6+6 posizioni orizzontli.

HORIZONTAL TEST TUBE RACKS AND THERMOMETER HOLDER

Made of polypropylene - Autoclavable. 6 + 6 horizontal positions.